



TZ MEDICAL PANNUS RETENTION SYSTEM

Instruction for Use

CZ-PRS-02

2 - 5.5" x 7.5" Pads

1 - 8' Straps



- 1) Peel backing from first adhesive pad and apply to underside of patient panniculus (lateral to belly button) on side where femoral access is desired. Ensure pad fully adheres flat against skin and that there are no bubbles or creases in pad. Position pad according to physician preference.
- 2) Peel backing from second adhesive pad and place pad on the underside/edge at the head of procedure table on opposite side of patient pad (at underside corner of table).
- 3) Place one end of strap (soft side down) on hook and loop surface of pad on patient panniculus, and drape remainder of strap over patient shoulder toward pad on corner head of procedure table.
- 4) Use one hand to grasp end of strap and pull it up toward pad on table while using the other hand to gently lift panniculus up into position.
- 5) Once desired positioning of patient panniculus has been achieved, attach remaining end of strap to hook and loop surface of pad on head of procedure table to secure.
- 6) To remove PRS from patient, remove strap from hook & loop pads and peel pad off of procedure table. To remove hook & loop pad from patient, start at one corner and gently peel edge of pad straight up with one hand, while holding patient skin down with the other hand. Discard PRS according to hospital protocol.

Step 1



Step 2



Step 3



Step 4



Step 5



* PRS pads may be left on patient up to 24 hours, and can be used to retract panniculus at conclusion of procedure to ensure surgical site remains dry and clean.

** If re-positioning the panniculus becomes necessary during the procedure, circulator/anesthesiologist will remove strap from hook & loop pad at head of procedure table. Adjust patient and panniculus positioning as necessary, then re-affix strap to hook & loop pad on underside of procedure table to re-secure.