


MicroLite8™ Radiation Protective Eyewear

Revolutionary 8-Base Lens Curve for Maximum Field of Vision and Lateral Protection in High Risk Surgical Environments



MicroLite8 features an exclusive 8-base lens curve for maximum visual detail and clarity at all angles; plus a contoured frame architecture to enhance lateral protection from scatter radiation and airborne contaminants.



Industry-Leading Features

- 8-Base lens curve for maximum field of vision and clarity.
- MicroLite Frame Architecture: comfortable facial contour design with integrated Optibloc™ side shields for maximum protection.
- SoftGrip™ nosepiece for comfortable fit.
- Weighs just 2.63 ounces, with 95% radiation reduction at 125 KVP Direct Beam.
- Pb Equivalent Lens: 0.75 mm (+/- .05) Pb glass.
- ProLiteHT™ high light transmission leaded lenses for unsurpassed optical clarity.
- ProLiteHT lenses are the lightest available anywhere, for maximum comfort during lengthy procedures.
- Curved Uniflex temples for flexible, contoured fit and enhanced protection.
- Soft, thermopolymer temple inserts for a comfortable grip.
- Widest range of prescriptions available, with fastest Rx turnaround in the industry
- One year warranty from date of sale.
- Accessory bonus: protective case, microfiber bag, anti fog lens cleaning solution, adjustable ClickFit retaining cord.

MicroLite8 midnight blue, Model FB108B (above) and silver, Model FB108G (right).



Exclusive!

Only MicroLite® Eyewear has ProLiteHT™ Lenses

- Up to 20% more light transmission than all other leaded eyewear brands.
- Provides ultimate contrast and visual detail in the operating room.

AttenuTech is the only company that designs and manufactures eyewear specifically for your protection in fluoroscopic environments. All MicroLite® Series radiation protective eyewear is prescriptionable by design. Our in-house lab will process and ship your single vision or bifocal prescription faster than anyone in the industry.

Your Authorized Dealer:



A Family of Medical Products Sparked by Your Ideas

7272 SW Durham Road, Ste #800
Portland, OR 97224
Tel: 800.944.0187 Fax: 503.639.0239
www.tzmedical.com