

MicroLiteM3™

Radiation Protective Eyewear

Revolutionary Adjustable Design for a Custom Fit

- **Sleek frame architecture** incorporates M3 Alloy, a proprietary blend of three metal compounds for lightweight durability
- **Spring hinge frame** for secure, precise fit
- **Adjustable SoftGrip™ nosepiece** conforms to your features for maximum comfort
- **MicroLiteM3** weighs just 2.2 oz.
- **Rapid Prescription Fulfillment:** The fastest turnaround in the industry



MicroLiteM3
Model FB103
Frame color: Grey

Additional MicroLiteM3 Features

- **ProLiteHT™** high light transmission leaded lenses for unsurpassed optical clarity
- **Integrated Optibloc™** leaded side shields provide lateral protection against scatter radiation
- **Attenuation:** 95% radiation reduction at 125 KVP Direct Beam
- **Pb Equivalent Lens:** 0.75 mm (+/- .05) Pb glass.
- **Widest range** of prescriptions available
- **One year warranty** from date of sale
- **Accessory bonus:** protective case, microfiber bag, anti-fog lens cleaning solution, adjustable ClickFit™ retaining cord.

Spring hinge frame and adjustable nosepiece put you in control of the articulation, affording a perfect three-point fit:

Protection • Lens Clarity • Comfort

AttenuTech is the only company that designs and manufactures eyewear specifically for your protection in fluoroscopic environments. All MicroLite® Series radiation protective eyewear is prescriptionable by design. Our in-house lab will process and ship your single vision or bifocal prescription faster than anyone in the industry.

Your Authorized Dealer:



A Family of Medical Products Sparked by Your Ideas

7272 SW Durham Road, Ste #800
Portland, OR 97224

Tel: 800.944.0187 Fax: 503.639.0239

www.tzmedical.com

Exclusive! **Only MicroLite® Eyewear** **has ProLiteHT™** **Lenses**

- Up to 20% more light transmission than all other leaded eyewear brands
- Provides ultimate contrast and visual detail in the operating room
- The lightest weight lenses available anywhere