

February 15, 2006

To Whom It May Concern:

TZ Medical certifies that the Lifelinks defibrillation electrodes have been validated with and considered compatible with the following defibrillators:

**Lifelinks Product Number:** P-111-T1, P-211-T1, P-121-T1, P-214-T1, P-122-T1


**Defibrillators:** Hewlett Packard Codemaster, HP 76000, HP 78000, Zoll PD 1200, Zoll PD 1400, Zoll M Series, Siemens 410, PhysioControl LP 5, PhysioControl LP 6, PhysioControl LP 8, PhysioControl LP 9, PhysioControl LP9P, PhysioControl LP 10, and PhysioControl 12

Currently, TZ Medical is in the process of verifying the compatibility of the TZ Medical defibrillation electrodes with the PhysioControl LP 20.

This validation documentation was submitted with our 510(K) number: **983196**

Additional validation documentation is located in the document control center at TZ Medical headquarters.

Sincerely,



Valerie Meyer  
TZ Medical Inc.