





# TZ Medical TruBrite

## Adhesive Retractor Light



### Product Features:

- Sterile, compact & disposable design that weighs less than 4 grams
- Provides continuous light for over 2 hrs
- Adjustable LED light, providing optimal visibility for any procedure
- Self-adhesive for effortless mounting on multiple surfaces
- Ideal for pacemaker pockets with Army-Navy retractors



### Product Info

<b>TruBrite</b>		
Item #	Description	UOM
ARL-1010	Adhesive Retractor Light	10/box