



# TZ MEDICAL COMFORT HALO

*Comfort Halo shown with:  
Comfort Cradle, CZ-100;  
Comfort Shield, CZ-005; &  
Patient Immobilizer, CZ-106*



*Comfort Halo  
Part# CZ-009*

## Comfort Halo

The TZ Medical Comfort Halo is designed to enhance safety and comfort for patients undergoing conscious sedation.

The board slides underneath the imaging table mattress, enhancing full view of the patient, and keeps the anesthesia drape off of the patient's face, improving patient satisfaction. The steel and vinyl tubing is flexible and can be molded into unlimited draping positions. The polycarbonate base is compatible with all bed types.

The Comfort Halo is ideal for all departments where radiographic imaging is used and sterile technique is essential.

## Product Features

- Increases patient comfort by keeping anesthesia drape off of patient's face during conscious sedation
- Offers an unobstructed view of patient during procedure
- Flexible design, may be molded to desired shape or height
- Visual reminder to perform "Time Out", enhancing patient safety/ improving staff compliance with JCAHO requirements
- Compatible with any bed type, no arm rails necessary

## For All Procedures Using Conscious Sedation Including:

- PPM/ICD/Bi-V ICD implants
- Chemo/Port placements
- CVP and PICC line insertion
- Antegrade femoral angiography
- Rt Heart Caths from Subclavian Approach



TZ Medical Inc.  
17750 SW Upper Boones Ferry Rd. #150  
Portland, OR 97224  
Office: 800.944.0187 Fax: 503.639.0239  
www.tzmedical.com

**TZ Medical**  
Sparked by your ideas

