



# TZ MEDICAL FLAT RADIAL BOARD



The dimensions of the Flat Radial Board are 44.4" long and 22.9" at its widest point, offering plenty of space and positioning options.

## TZ Medical Flat Radial Board

The TZ Medical Flat Radial Board is a lightweight carbon fiber board specifically shaped to allow convenient radial access in conjunction with medical imaging and patient comfort. By combining radial access and a landing zone for equipment, the Flat Radial Board takes radial access to a whole new level.

The Flat Radial Board can be used with the new Radial Runway, which helps keep the patient's wrist in position for radial access. Both the Flat Radial Board and the Radial Runway can be used for radial access on either the right or left side.

*Flat Radial Board shown in use with Radial Runway Part# RRU-008 & RRU-009*



Flat Radial Board  
Part# CZ-400

### Product Features

- Can be used on either the right or left side for radial access
- Large landing area for equipment
- Lightweight construction
- Invisible under fluoroscopy
- Can be used for all diagnostic and interventional imaging procedures

TZ Medical Inc.  
17750 SW Upper Boones Ferry Rd. #150  
Portland, OR 97224  
Office: 800.944.0187 Fax: 503.639.0239  
www.tzmedical.com

**TZ Medical**  
Sparked by your ideas

