



# TZ MEDICAL EZ HOLD & NEPTUNE DISCS

## TZ Medical EZ Hold

The TZ Medical EZ Hold is a manual-assist compression device developed to aide the clinician in manual femoral pressure.

The EZ Hold was designed and patented by an RN who sought to reduce hand fatigue and direct contact with bodily fluids. The ergonomic design of the EZ Hold permits the clinician to hold pressure comfortably with one hand while palpating the pulse with the other.

The TZ Medical EZ Hold is best used in conjunction with our line of Neptune Discs. These come in two styles, and have the advantage of having a Neptune pad built into the disc, which speeds up hemostasis and reduces holding times.



*Oval Universal Disc  
with Calcium Alginate  
Part# 8880-07*



*Round Universal Disc  
with Calcium Alginate  
Part# 8880-08*



*2 lb. EZ Hold  
Part# 8260-20*

## Neptune Discs

Femoral pressure discs ensure that just the right amount of pressure is applied to the access site when used in conjunction with either our EZ hold or QuicKlamp products.

TZ Medical Neptune Discs use the power of calcium alginate to speed up hemostasis. The calcium alginate is sustainably harvested from seaweed and is an all natural product.

TZ Medical Inc.  
17750 SW Upper Boones Ferry Rd. #150  
Portland, OR 97224  
Office: 800.944.0187 Fax: 503.639.0239  
www.tzmedical.com

**TZ Medical**  
Sparked by your ideas

