




The Halo Flipper is a versatile radiation shield designed to conform and secure to the TZ Medical Comfort Halo to enable increased radiation protection at the head of the procedure table. The Halo Flipper attenuates up to 83% of radiation at 90 kVp with a lead equivalence of 0.125mm. The Halo Flipper is designed with adhesive strips on the back side of the shield and can also be used as an additional radiation protection drape when adhered to the edge of the procedure table.



Features

- Lead-free and attenuates up to 83% at 90 kVp
- Provides additional radiation protection for staff members at the head of the procedure table
- Adhesive strips are designed to adhere securely around the flexible tubing of the Comfort Halo
- Can be used as an additional radiation protection table drape when adhered to the side of the procedure table
- Single-use and packaged sterile



Item Number	Description	UOM
RAD-1020	Halo Flipper Radiation Shield	Box of 5

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

23675 R0.1