


The Peripheral Flipper is a versatile radiation shield that adheres to the procedure drape during peripheral procedures for additional protection against scatter radiation. The device is designed with adhesive strips on the front and back, allowing users to place the shield in multiple locations for optimal protection. The Peripheral Flipper attenuates up to 83% of radiation at 90 kVp with a lead equivalence of 0.125mm.

## Features

- Lead-free and attenuates up to 83% at 90 kVp
- Adhesive strips on both sides help prevent slipping and allow the shield to be used in multiple locations
- Provides additional radiation protection during peripheral procedures
- Single-use and packaged sterile



Item Number	Description	UOM
RAD-1030	Peripheral Flipper Radiation Shield	Box of 5

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)  
23676 R0.1