




The TZ Medical Labor and Delivery Drape Holder is an adaptable solution that slides under any mattress and provides simple drape positioning for the comfort of both the patient and medical staff. Crafted with a durable polypropylene base and two adjustable 30" flexible goosenecks, the Drape Holder offers versatility in shaping to optimally position your patients and ensure their comfort.



## Features

- Eliminates the need for using bulky, large-footprint IV poles for drape positioning
- Reduces equipment congestion at the head of the bed and creates more room for medical staff and patient family
- Easy access for skin-to-skin with newborn
- Flexible gooseneck arms can be molded to any desired shape or height to accommodate patients of all shapes and sizes
- Compatible with all bed types; no arm rails or additional equipment necessary
- Sturdy, easy-to-clean polypropylene base
- Easily repositioned to adapt to changing procedural needs



Item Number	Description	UOM
SDH-0005	L&D Drape Holder	1 Each

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)