



TZ Medical now offers a line of stainless steel instrument tables. Each lightweight stainless steel table has heavy-duty swivel casters, a lower shelf, fully welded legs, and 16-gauge steel. Available in three lengths, TZ tables provide a stable, easy-to-clean surface for all procedures throughout the hospital.

Features

- 16-gauge stainless steel
- Closed hems for simple cleaning
- Convenient lower shelf for storage
- Comes in three lengths: 48", 60", and 72"
- Tables are delivered fully assembled
- Made in the USA



Item Number	Description	UOM
SSIT-48	Instrument Table - 48"L x 24"W x 34"H	1 Each
SSIT-60	Instrument Table - 60"L x 24"W x 34"H	1 Each
SSIT-72	Instrument Table - 72"L x 24"W x 34"H	1 Each

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

23982 R0.2