



TZ Medical
Sparked by your ideas

STARTROL MEDICAL LIGHTS



StarTrol LED Medical Lights produce a bright white light, (94 CRI and 4700K). They are ideal for true color rendition, and are cool and comfortable to work under, with virtually no radiant heat being generated by the LED's. The lights are energy-efficient, using less than one-third the energy of a halogen light, and require minimal maintenance. StarTrol LED Medical Lights are durable, precision machined in the USA and feature a 10-year warranty. TZ Medical offers a wide variety of light configurations, including ceiling-mounted, wall-mounted, and mobile rolling units, which are sure to fit into any workspace.

Features

- Energy Efficient & Eco-friendly
- Bright, natural white light
- 10-year limited warranty
- Cool to work under
- 40,000 hours of LED life
- Infinite variable positioning
- Multi-level dimming



Item Number	Description	UOM
ST-4X4-GSM	4 Light Spring Arm Mobile	Each
ST-8X4-GSM	8 Light Spring Arm Mobile	Each
ST-4X4-GSW	4 Light Spring Arm Wall	Each
ST-8X4-GSW	8 Light Spring Arm Wall	Each
ST-4X4-GS1	4 Light Spring Arm Single	Each
ST-8X4-GS1	8 Light Spring Arm Single	Each
ST-4X4-GS2	4 Light Spring Arm Dual	Each
ST-8X4-GS2	8 Light Spring Arm Dual	Each

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

15763 R3.0

Ceiling Mounted Units



TZ Medical
Sparked by your ideas

Dual Arm Ceiling Mounted Units



Wall Mounted Units



Mobile Units



FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

15763 R3.0